



**WATER
SAMPLING KIT**

WATER SAMPLING KIT

LABPLAS sterile **WATER SAMPLING KIT** is designed for collecting, transporting and analyzing chlorinated and non-chlorinated water samples in aseptic conditions. It contains a non-toxic and neutralizing tablet of active sodium thiosulfate, to allow detection of potential pathogens. UPA compliant for the analysis of drinking water, it is an economical and efficient alternative to rigid sampling containers.

A quick and convenient way to collect water samples!

Features and benefits

- Made of FDA compliant virgin polyethylene
- Graduated lines on the bags for different volumes
- Can be transported in a cooler
- Adequate for sampling controls of the FDA, USDA, HACCP and CFIA
- Sterilization process done by gamma irradiation




Inspired by your challenges

LABPLAS

VISIT LABPLAS.COM FOR ADDITIONAL INFORMATION ABOUT OUR PRODUCTS | Tel: 1.450.649.7343



WATER SAMPLING KIT BAG INSTRUCTION



Inspired by your challenges

LABPLAS

VISIT LABPLAS.COM FOR ADDITIONAL INFORMATION ABOUT OUR PRODUCTS | Tel: 1.450.649.7343

