



Sani-S'wipe

SANI-S'WIPE

The SANI-S'WIPE is designed to detect the presence of biological contamination such as *Listeria*, *Salmonella*, *E. coli* and other food-borne pathogens that may be found on various surfaces in contact with contaminated materials or on the product itself. The sampling kit is ideal for the control of large areas (over 100 cm²).

Easy to use, Sani-S'wipe relies on a quick and effective gesture to ensure that your controlled surface is representative!

Features and benefits

- The kit complies with ISO 18593
- Ensures a representative sample
- Reduced sample collection time compared to excision methods
- Laser scoring opening system by laser marking line to guarantee the seal
- The kit can be used in the homogenizer
- Available with a dry or pre-moisturized cloth with a neutralizing or non-neutralizing solution
- High-quality cloth that does not discolour during sampling and analysis




Inspired by your challenges

LABPLAS

VISIT LABPLAS.COM FOR ADDITIONAL INFORMATION ABOUT OUR PRODUCTS | Tel: 1.450.649.734



SANI-S'WIPE BAG INSTRUCTION



Inspired by your challenges

LABPLAS

VISIT LABPLAS.COM FOR ADDITIONAL INFORMATION ABOUT OUR PRODUCTS | Tel: 1.450.649.734

