

## GENERAL INFORMATION

Category : Sterile sampling kits

Family : Sani-Sponge

Lifespan: 5 years

## TECHNICAL DESCRIPTION

This kit consists of a Twirl'em bag and a biocide-free sponge for sampling various surfaces.



## Volume guide

<sup>1</sup>Total volume: The total volume corresponds to the maximum capacity the bag can hold when it is filled to the brim, without the possibility of closing it.

<sup>2</sup>Closure volume: The closure volume corresponds to the maximum capacity the bag can hold when it is closed with a triple twist to ensure safe transport.

<sup>3</sup>Functional volume: The functional volume corresponds to the volume of diluent or sample according to the test requirements.

## SPECIFIC INFORMATION

## CONTAINER

Type of container :	Bag
Material :	Polyethylene blend
Color :	Transparent
Dimension :	114 x 229 mm / 4.5 x 9 in
Thickness :	2.5 mil.in/ 63.5 micron / 0.0635 mm
Total volume <sup>1</sup> :	650 ml / 22 oz
Closure volume <sup>2</sup> :	390 ml / 13 oz
Functional volume <sup>3</sup> :	100 ml / 3 oz
Printing type :	Writing area
Opening system :	Perforated line
Closing system :	Attachment with 2 round wires
Sterile :	Yes
End of product life :	Recyclable

## INSERTION

Insertion type :	Cellulose sponge
Description :	Cellulose sponge. Yellow color. Dimensions (L x W): 3" x 1.5" (77 mm x 39 mm). Free from biocides.
Biodegradable :	Non-recyclable

## INSERTION

Insertion type :	
Description :	
Biodegradable :	

## SOLUTION

Name :	NA
Type of solution :	NA
pH :	NA

## OTHER

--

## AVAILABLE DOCUMENTS

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

## PACKAGING INFORMATION

Outer box dimension : (W x D x H)	9.50 X 14.13 X 16.25
Box weight :	7.80 LB / 3.54 KG
Conditioning :	400 (4 x100)
Storage condition :	Store in a dry place at room temperature

Reach out to us for additional resources, if applicable to this product.

## DECLARATION

CFIA	LABPLAS sampling bags are a solution that may be used in the CFIA Preventive Control Plan (PCP) for the seven principles of the HACCP system. The PCP is a Canadian federal initiative, under the Safe Food for Canadians Regulations (SFCR).
EU	The materials used to manufacture LABPLAS sampling bags meet, where applicable, the Eu No10/2011 standards for food contact with respect to particle migration.
FDA	The plastic film used in the manufacture of the LABPLAS sampling bag meets the requirements of 21 CFR 177.1520 of the Food and Drug Administration.
Biocide-free	Biocide-free material
Sterile	Sterility is provided by irradiation. The sterilization dose ensuring a sterility assurance level (SAL) of 10 <sup>-3</sup> has been established according to ISO 11137.

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com