

GENERAL INFORMATION

Category : Sterile sampling kits

Family : Wipa-Boots

Lifespan: 24 months

TECHNICAL DESCRIPTION

A kit consisting of sterile woven cotton slippers inserted in a sterile bag. The kit allows for large-area sampling at ground level.



Volume guide

¹Total volume: The total volume corresponds to the maximum capacity the bag can hold when it is filled to the brim, without the possibility of closing it.

²Closure volume: The closure volume corresponds to the maximum capacity the bag can hold when it is closed with a triple twist to ensure safe transport.

³Functional volume: The functional volume corresponds to the volume of diluent or sample according to the test requirements.

SPECIFIC INFORMATION

CONTAINER

Type of container :	Bag
Material :	PE linear tear film
Color :	Transparent
Dimension :	160 x 279 mm / 6.3 x 11 in
Thickness :	4 mil.in / 101.6 micron / 0.1016 mm
Total volume ¹ :	1200 ml / 41 oz
Closure volume ² :	740 ml / 25 oz
Functional volume ³ :	250 ml / 8 oz
Printing type :	Form
Opening system :	Notch
Closing system :	Attachment with 2 round wires
Sterile :	Yes
End of product life :	Recyclable

INSERTION

Insertion type :	Pedichiffonnette
Description :	Cotton knitted towel. Size: one size. Free from biocides
Biodegradable :	Non-recyclable

INSERTION

Insertion type :	
Description :	
Biodegradable :	

SOLUTION

Name :	Tryptone-Salt Broth
Type of solution :	Transport medium, non-neutralizing
pH :	7.00±0.2

OTHER

--	--

AVAILABLE DOCUMENTS

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

PACKAGING INFORMATION

Outer box dimension : (W x D x H)	15.25 X 22.75 X 7.00
Box weight :	12.26 LB / 5.56 KG
Conditioning :	100 (1 x 100)
Storage condition :	Upon receipt, store between 2 °C and 8 °C

Reach out to us for additional resources, if applicable to this product.

DECLARATION

CFIA	LABPLAS sampling bags are a solution that may be used in the CFIA Preventive Control Plan (PCP) for the seven principles of the HACCP system. The PCP is a Canadian federal initiative, under the Safe Food for Canadians Regulations (SFCR).
EU	The materials used to manufacture LABPLAS sampling bags meet, where applicable, the Eu No10/2011 standards for food contact with respect to particle migration.
FDA	The plastic film used in the manufacture of the LABPLAS sampling bag meets the requirements of 21 CFR 177.1520 of the Food and Drug Administration.
Biocide-free material	Biocide-free material for detection of salmonella
Sterile	Sterility is provided by irradiation. The sterilization dose ensuring a sterility assurance level (SAL) of 10 ⁻³ has been established according to ISO 11137.

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com