

## GENERAL INFORMATION

Category : Sterile sampling kits

Family : Sani-Stick

Lifespan: 15 months

## TECHNICAL DESCRIPTION

This kit consists of a Twirl'em bag and a biocide-free sponge on a handle designed for vigorous sampling of various surfaces and reaching difficult areas.



## SPECIFIC INFORMATION

## CONTAINER

Type of container :	Bag
Material :	Thermoplastic laminate
Color :	Transparent
Dimension :	190 x 292mm / 7.5 X 11.5 po
Thickness :	3.5 mil.in / 0.0889 mm / 88.9 micron
Total volume :	2100 ml / 70 oz
Functional volume :	1200 ml / 41 oz
Printing type :	Writing area
Opening system :	Laser scoring line
Closing system :	Attachment with 2 round wires
Sterile :	Yes
End of product life :	Garbage

## INSERTION

Insertion type :	Polyurethane sponge
Description :	Hydrophilic polyurethane foam sponge. White colour. Dimensions (L x W): 3" x 1.25" (77 mm X 31.8 mm). Free from biocides.
Biodegradable :	Garbage

## INSERTION

Insertion type :	Rigid Handle
Description :	Rigid polypropylene handle assembled with sponge. Red color. Equipped with a contactless sponge release mechanism
Biodegradable :	Recyclable

## SOLUTION

Name :	Dey-Engley neutralizing broth
Type of solution :	Neutralizing
pH :	7.60±0.2

## OTHER

Sterile blue gloves are the perfect tool to complement your sampling bags and kits. The Blue polyethylene gloves are packed in individual pairs, and then irradiated to sterilize them. They are folded and packed optimally to maintain sterility.

**AVAILABLE DOCUMENTS**

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

**PACKAGING INFORMATION**

Outer box dimension : (W x D x H)	13.63 X 17.25 X 14.00
Box weight :	8.50 LB / 3.86 KG
Conditioning :	100 (10 x 10)
Storage condition :	Upon receipt, store between 2 °C and 8 °C

Reach out to us for additional resources, if applicable to this product.

**DECLARATION**

CFIA	LABPLAS sampling bags are a solution that may be used in the CFIA Preventive Control Plan (PCP) for the seven principles of the HACCP system. The PCP is a Canadian federal initiative, under the Safe Food for Canadians Regulations (SFCR).
EU	The materials used to manufacture LABPLAS sampling bags meet, where applicable, the Eu No10/2011 standards for food contact with respect to particle migration.
FDA	The plastic film used in the manufacture of LABPLAS sampling bags complies with 21 CFR 177.1520 [(c) 3.2c] of the Food and Drug Administration. Its use is limited to temperatures below 212 F according to Table 2 of 21 CFR 176.170 (c).
Biocide-free	Biocide-free material
Sterile	Sterility is provided by irradiation. The sterilization dose ensuring a sterility assurance level (SAL) of 10 <sup>-3</sup> has been established according to ISO 11137.

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com