

## GENERAL INFORMATION

Category: Racks and complementary tools

Family: Work station

Lifespan: 27 years

## TECHNICAL DESCRIPTION

The Work station keeps the sampling bag open in the vertical position.



## SPECIFIC INFORMATION

ITEM	Rack
Material :	Polycarbonate
Color :	Transparent
Resistance :	Temperature resistance from -100°C to 135°C
Dimension :	127 x 133 x 270 mm / 5 x 5.25 x 10.5 in
Autoclavable :	No
Sterile :	No
End of product life :	Recyclable

## PACKAGING INFORMATION

Outer box dimension : (W x D x H)	8.06 X 13.13 X 6.63
Box weight:	1.00 LB / 0.45 KG
Conditioning:	1 (1x1)
Storage condition:	No storage conditions

## OTHER

**AVAILABLE DOCUMENTS**

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

Reach out to us for additional resources, if applicable to this product.

**DECLARATION**

Non-sterile	This product is offered non-sterile.
-------------	--------------------------------------

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com