

GENERAL INFORMATION

Category: Racks and complementary tools

Family: Dispenser

Lifespan: 27 years

TECHNICAL DESCRIPTION

The dispenser is a tool that allows the distribution of a roll of 250 bags, it is resistant and compact.



SPECIFIC INFORMATION

ITEM	Tool
Material :	Stainless steel
Color :	Steel gray
Resistance :	Temperature resistance -80°C to 80°C
Dimension :	203 x 114 x157 mm / 8 x 4.5 x 6 in
Autoclavable :	Yes
Sterile :	No
End of product life :	Recyclable

PACKAGING INFORMATION

Outer box dimension : (W x D x H)	6.75 X 8.00 X 7.75
Box weight:	0.30 LB / 0.14 KG
Conditioning:	1 (1x1)
Storage condition:	No storage conditions

OTHER

AVAILABLE DOCUMENTS

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

Reach out to us for additional resources, if applicable to this product.

DECLARATION

Non-sterile	This product is offered non-sterile.
-------------	--------------------------------------

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com