

GENERAL INFORMATION


Category: Racks and complementary tools

Family: LABrack

Lifespan: 27 years

TECHNICAL DESCRIPTION

The LABrack is a stainless steel rack that allows you to carry multiple bags easily and safely. It can go in the incubator and autoclave, and can be stacked for space saving.



SPECIFIC INFORMATION

ITEM	Rack
Material :	Stainless steel
Color :	Steel gray
Resistance :	Temperature resistance -80°C to 80°C
Dimension :	133 mm / 5.25 in
Autoclavable :	Yes
Sterile :	No
End of product life :	Recyclable

PACKAGING INFORMATION

Outer box dimension : (W x D x H)	9.63 X 16.38 X 9.63
Box weight:	7.00 LB / 3.18 KG
Conditioning:	1 (1x1)
Storage condition:	No storage conditions

OTHER

This product comes with 10 closure clips, they are used to close bags that do not have a closure system, The closure clip is made up of a cylinder and a piece in the shape of a half-moon that one clips/closes by pressure on each other.

AVAILABLE DOCUMENTS

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

Reach out to us for additional resources, if applicable to this product.

DECLARATION

Non-sterile	This product is offered non-sterile.
-------------	--------------------------------------

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com