

## GENERAL INFORMATION


Category: Racks and complementary tools

Family: LABrack

Lifespan: 27 years

## TECHNICAL DESCRIPTION

The LABrack is a stainless steel rack that allows you to carry multiple bags easily and safely. It can go in the incubator and autoclave, and can be stacked for space saving.



## SPECIFIC INFORMATION

ITEM	Rack
Material :	Stainless steel
Color :	Steel gray
Resistance :	Temperature resistance -80°C to 80°C
Dimension :	197 x 419 mm / 7.8 x 16.5 in
Autoclavable :	Yes
Sterile :	No
End of product life :	Recyclable

## PACKAGING INFORMATION

Outer box dimension : (W x D x H)	9.63 X 16.38 X 9.63
Box weight:	7.00 LB / 3.18 KG
Conditioning:	1 (1 x 1)
Storage condition:	No storage conditions

## OTHER

This product comes with 10 closure clips, they are used to close bags that do not have a closure system, The closure clip is made up of a cylinder and a piece in the shape of a half-moon that one clips/closes by pressure on each other.

## AVAILABLE DOCUMENTS

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

Reach out to us for additional resources, if applicable to this product.

## DECLARATION

Non-sterile	This product is offered non-sterile.
-------------	--------------------------------------

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com