

## GENERAL INFORMATION

Category: Racks and complementary tools

Family: Rack'N'Lab

Lifespan: 27 years

## TECHNICAL DESCRIPTION

The Rack'N'Lab is a stand designed to facilitate the transport of bags and to hold filled bags vertically. It can go in the cooler, incubator and autoclave, and can be stored flat quickly.



## SPECIFIC INFORMATION

ITEM	Rack
Material :	Polypropylene
Color :	White
Resistance :	Temperature resistance from 0°C to 135°C
Dimension :	305 x 171 x 76 mm / 12 x 6.75 x 3 in
Autoclavable :	Yes
Sterile :	No
End of product life :	Recyclable

## PACKAGING INFORMATION

Outer box dimension : (W x D x H)	1.25 X 13.00 X 23.00
Box weight:	2.20 LB / 1.00 KG
Conditioning:	3 (1 x 3)
Storage condition:	No storage conditions

## OTHER

**AVAILABLE DOCUMENTS**

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

Reach out to us for additional resources, if applicable to this product.

**DECLARATION**

Non-sterile	This product is offered non-sterile.
-------------	--------------------------------------

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com