

GENERAL INFORMATION

Category: Racks and complementary tools

Family: Rack'N'Lab

Lifespan: 27 years

TECHNICAL DESCRIPTION

The Rack'N'Lab is a stand designed to facilitate the transport of bags and to hold filled bags vertically. It can go in the cooler, incubator and autoclave, and can be stored flat quickly.



SPECIFIC INFORMATION

ITEM	Rack
Material :	Polypropylene
Color :	White
Resistance :	Temperature resistance from 0°C to 135°C
Dimension :	305 x 171 mm / 12 x 6.8 in
Autoclavable :	Yes
Sterile :	No
End of product life :	Recyclable

PACKAGING INFORMATION

Outer box dimension : (W x D x H)	1.25 X 13.00 X 23.00
Box weight:	2.20 LB / 1.00 KG
Conditioning:	3 (1 x 3)
Storage condition:	No storage conditions

OTHER

AVAILABLE DOCUMENTS

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

Reach out to us for additional resources, if applicable to this product.

DECLARATION

Non-sterile	This product is offered non-sterile.
-------------	--------------------------------------

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com