

## GENERAL INFORMATION

Category: Racks and complementary tools

Family: Sterile blue gloves

Lifespan: 27 years

## TECHNICAL DESCRIPTION

Made of polyethylene, the Sterile blue gloves are packaged in pairs and folded so that they can be slipped on safely.



## SPECIFIC INFORMATION

ITEM	Tool
Material :	Polyethylene
Color :	Blue
Resistance :	Temperature resistance -80°C to 80°C
Dimension :	0 x 0 mm / 0 x 0 in
Autoclavable :	No
Sterile :	Yes
End of product life :	Recyclable

## PACKAGING INFORMATION

Outer box dimension : (W x D x H)	7.88 X 8.13 X 1.63
Box weight:	1.50 LB / 0.68 KG
Conditioning:	50 (1 x 50p)
Storage condition:	No storage conditions

## OTHER

## AVAILABLE DOCUMENTS

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

Reach out to us for additional resources, if applicable to this product.

## DECLARATION

Sterile	Sterility is provided by irradiation. The sterilization dose ensuring a sterility assurance level (SAL) of $10^{-3}$ has been established according to ISO 11137.
---------	---

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com