



## GENERAL INFORMATION

Category: Racks and complementary tools

Family: Reach'em

Lifespan: 27 years

## TECHNICAL DESCRIPTION

The Reach'em is a lightweight yet durable adapter designed to help the Sani-Stick reach hard-to-access surfaces. It attaches securely to a standard telescopic pole, extending reach during sampling operations.



## SPECIFIC INFORMATION

ITEM	Tool
Material :	Polylactic Acid (PLA)
Color :	Green
Resistance :	Temperature resistance -80°C to 80°C
Dimension :	36 x 118 mm / 1.4 x 4.6 in
Autoclavable :	No
Sterile :	No
End of product life :	Recyclable

## PACKAGING INFORMATION

Outer box dimension : (W x D x H)	5.75 X 8.88 X 2.63
Box weight:	0.23 LB / 0.10 KG
Conditioning:	1 (1 x 1)
Storage condition:	No storage conditions

## OTHER

## AVAILABLE DOCUMENTS

Data Sheet	Certificate of Compliance
Certificate of Analysis	Safety Data Sheet (SDS)
Certificate of Sterility	Pyrogen Declaration
DNase/RNase	

Reach out to us for additional resources, if applicable to this product.

## DECLARATION

Non-sterile	This product is offered non-sterile.
-------------	--------------------------------------

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com