

## GENERAL INFORMATION

Category: Racks and complementary tools

Family: Reach'em

Lifespan: 27 years

## TECHNICAL DESCRIPTION

The Reach'em is a lightweight yet durable adapter designed to help the Sani-Stick reach hard-to-access surfaces. It attaches securely to a standard telescopic pole, extending reach during sampling operations.



## SPECIFIC INFORMATION

| ITEM                  | Tool                                 |
|-----------------------|--------------------------------------|
| Material :            | Polylactic Acid (PLA)                |
| Color :               | Green                                |
| Resistance :          | Temperature resistance -80°C to 80°C |
| Dimension :           | 36 x 118 mm / 1.4 x 4.6 in           |
| Autoclavable :        | No                                   |
| Sterile :             | No                                   |
| End of product life : | Recyclable                           |

## PACKAGING INFORMATION

|                                      |                       |
|--------------------------------------|-----------------------|
| Outer box dimension :<br>(W x D x H) | 7.50 X 11.00 X 0.00   |
| Box weight:                          | 0.23 LB / 0.10 KG     |
| Conditioning:                        | 1 (1 x 1)             |
| Storage condition:                   | No storage conditions |

## OTHER

## AVAILABLE DOCUMENTS

|                          |                           |
|--------------------------|---------------------------|
| Data Sheet               | Certificate of Compliance |
| Certificate of Analysis  | Safety Data Sheet (SDS)   |
| Certificate of Sterility | Pyrogen Declaration       |
| DNase/RNase              |                           |

Reach out to us for additional resources, if applicable to this product.

## DECLARATION

|             |                                                                                 |
|-------------|---------------------------------------------------------------------------------|
|             | PLA material declared compliant with FDA 21 CFR 175.300 (indirect food contact) |
| Non-sterile | This product is offered non-sterile.                                            |
|             | Material compliant with Directive 2011/65/EU                                    |

1951, rue Nobel, Sainte-Julie, Québec, Canada, J3E 1Z6

Tél : 1.450.649.7343 | Fax : 1.450.649.3113

labplas@labplas.com | www.labplas.com